



# ***Adaptive Coatings in Micro-Capillaries and Micro-Channels as New Means to Develop Integrated Micro-Fluidic Optical Sensors***

Larisa Florea

October 2011

# Outline

## I. Introduction

- optical sensors in microfluidics
- stimuli-responsive materials

## II. Polyaniline functionalised micro-fluidic channels for:

- pH sensing
- study mixing

## III. Spiropyran functionalised micro-capillaries for:

- solvent sensing




Outline

I. Introduction

II. Polyaniline functionalised  
micro-channels

III. Spiropyran functionalised  
micro-capillaries

# Optical Sensors in Microfluidics



I.

Optical Chemical Sensors

Stimuli responsive materials

# Stimuli-responsive Materials

Materials whose characteristics can be changed using an external stimulus

Thermal

Mechano

Photo

pH

Chemo

I.




Optical Chemical Sensors

Stimuli responsive materials



## II. Polyaniline functionalised micro-fluidic channels

- pH sensing
- study mixing



II.

Background

Polyaniline  
Micro-channels


pH Sensing

Mixing Process

Conclusions


# Polyaniline Nanofibres

- low cost, easy synthesis
- reversible acid-base doping-dedoping chemistry
- environmental stability



Insulating State  
Blue/Violet Colour

Conducting State  
Green



J.X. Huang, S. Viril, B.H. Weller, R.B. Kaner / J.Am.Chem.Soc. 125 (2003), 314-315

II.

Background

Polyaniline  
Micro-channels


pH Sensing

Mixing Process


Conclusions

# Micro-chip fabrication


✓PDMS is poured onto master mold, cured at 80°C for 1 h and removed from mold.



✓PDMS and glass slide are treated with oxygen plasma.



✓PDMS and glass slide are brought together.



L. Yu, C.M. Li, Y. Liu *et al.* / Lab Chip, 9 (2009), 1243–1247.

II.

Background



Polyaniline  
Micro-channels

pH Sensing

Mixing Process

Conclusions

# Micro-channel functionalisation process



II.

Background

Polyaniline  
Micro-channels


pH Sensing

Mixing Process


Conclusions

# Micro-channels


500 $\mu$ m x 1000 $\mu$ m



1000 $\mu$ m x 100 $\mu$ m



45 $\mu$ m x 50 $\mu$ m



II.

Background

Polyaniline  
Micro-channels


pH Sensing

Mixing Process

Conclusions



# Fast Response



II.

Background


Polyaniline  
Micro-channels

pH Sensing


Mixing Process

Conclusions


# pH sensing in continuous flow



pH 2




pH 4





pH 12

Dedoping process



Emeraldine Salt (ES)



Emeraldine Base (EB)

II.




Background

Polyaniline  
Micro-channels

pH Sensing

Mixing Process

Conclusions



L. Florea, E. Lahiff, D. Diamond, F. Benito-Lopez / MicroTas 2011, shortlisted for best poster.

L. Florea, E. Lahiff, D. Diamond, F. Benito-Lopez / Lab Chip, to be submitted.

**II.**

Background

Polyaniline  
Micro-channels

pH Sensing


Mixing Process

Conclusions



# Dynamic pH sensing

NaOH  $10^{-2}$  M



HCl  $10^{-2}$  M

II.

Background



Polyaniline  
Micro-channels

pH Sensing

Mixing Process

Conclusions

# Mixing process



L. Florea, E. Lahiff, D. Diamond, F. Benito-Lopez / MicroTas 2011, shortlisted for best poster.

**II.**

Background



Polyaniline  
Micro-channels

pH Sensing

Mixing Process

Conclusions

# Mixing process



II.

Background

Polyaniline  
Micro-channels

pH Sensing

Mixing Process

Conclusions

# Conclusions

- First example of polyaniline coated micro-channels
- Self-diagnostic for continuous flow device
- Simple and fast photometric method to measure pH
- Replace the glass layer with glass-ITO -> electro-chemical sensing of redox active species

II.

Background


Polyaniline  
Micro-channels

pH Sensing

Mixing Process

Conclusions

### III. Spiropyran functionalised micro-capillaries for solvent sensing



III.

Background


SP-monomer

Sp-polymer  
Brushes

Solvent Sensing


Conclusions

# Spiropyran



A : Spiropyran SP (closed, colorless)

B : Merocyanine MC (open, colored)



III.

Background


SP-monomer

Sp-polymer  
Brushes


Solvent Sensing

Conclusions





**Fig 1.** Experimental setup for measurement of capillary rise in a spiropyran monolayer functionalised capillary [1].



**Fig 2.** PDMS/glass hybrid micro-fluidic device functionalised with a monolayer of spiropyran molecules [2].

## Low loading of spiropyran molecule

[1] R. Rosario et al. / Langmuir, 18 (2002) 8062-8069

[2] F. Benito-Lopez et al. / Sensors and Actuators B 140 (2009) 295–303

III.


Background

SP-monomer

Sp-polymer  
Brushes

Solvent Sensing

Conclusions



↑ - spiropyran molecule

- polymer brushes
- high loading of spiropyran molecule
- 3D arrangement

Micro-capillary : Convenient platform for rapid analysis and detection

## Advantages

- act as a mechanical support for the optically sensitive layer
- represents an optical waveguide structure
- suitable for real-time continuous flow analysis
- requires very small volume of analyte

III.

Background



SP-monomer

Sp-polymer  
Brushes

Solvent Sensing

Conclusions





III.

Background

SP-monomer


Sp-polymer  
Brushes

Solvent Sensing


Conclusions

# Spiropyran polymeric brushes in micro-capillaries

## Silanisation



## Attachment of the catalyst



III.

Background


SP-monomer

Sp-polymer  
Brushes

Solvent Sensing

Conclusions

# Spiropyran polymeric brushes in micro-capillaries



III.

Background


SP-monomer

Sp-polymer  
Brushes

Solvent Sensing

Conclusions


# Characterisation



↓ UV light



↓ Remove mask



↓ White light



L. Florea, A. Hennart, D. Diamond, F. Benito-Lopez / Sens. Actuators B, submitted.

III.



Background

SP-monomer

Sp-polymer  
Brushes

Solvent Sensing

Conclusions



III.

Background

SP-monomer


Sp-polymer  
Brushes

Solvent Sensing

Conclusions

# Solvatocromic Properties

The colour of the MC form depends on the difference in polarity between the photo-excited MC form and the conjugated zwitterionic ground state



- The absorption band of MC form undergoes a hypsochromic (blue) shift in solvents of increasing polarity (negatively solvatochromism).

U.I. Minkin / Chem. Reviews, 104 (2004) 2751-2776.

III.

Background



SP-monomer

Sp-polymer  
Brushes

Solvent Sensing

Conclusions

# Solvent Detection



L. Florea, A. Hennart, D. Diamond, F. Benito-Lopez / Chem Comm., submitted.

III.

Background

SP-monomer

Sp-polymer  
Brushes

Solvent Sensing

Conclusions

# Conclusions

- First example of spiropyran polymer brushes in micro-capillaries
- Self-diagnostic for continuous flow device
- Solvent detection in micro-capillaries
- Sensing behaviour can be switched on/off remotely using light



# Acknowledgments

- Dr. Fernando Benito-Lopez
- Alexandre Hennart  
University of Namur, Belgium
- Dr. Kevin Fraser
- Dr. Emer Lahiff
- Prof. Dermot Diamond
- IRCSET – Embark Initiative
- CLARITY – award (07/CE/ I1147)



# Thank you for your attention !!

