

Real-time Analysis of Sweat Using Integrated Chemical Sensors

S. Coyle, V. Curto, F. Benito Lopez, D. Diamond

- Introduction
 - ✓ Sweat Analysis
 - ✓ Wearable sensors
- Sodium Sensors
- pH Sensors
- Lactate Sensor
- Conclusions

The sweat test is routinely used for the diagnosis of cystic fibrosis

Other conditions reported to cause a high salt level in sweat:

- Adrenal glands malfunction (adrenal insufficiency or Addison's disease)
- Hypothyroidism
- Diabetes Insipidus (caused by inability of the body to retain water)
- Kidney failure

Sweat as a removal pathway for drugs

- Legal (salicylic acid, antipyrine, phenol) and illegal (cocaine, heroine, amphetamines) drugs
- Forensic medicine is trying to exploit this to detect drug abuse

Clinical Medicine:

Chronic disease management (heart, renal failure)

Peripheral vascular disease (therapy efficacy)

Post-surgery management (hydration protocols)

Physical Training:


Personalised training protocols and reintegration

Down's syndrome children training program

People with thermoregulatory impairment

Technology comparison

Whole body sweat collection techniques



S. M. Shirreffs, *et al.*, J Appl Physiol 82 (1997) 336-34





Minimum discomfort





PharmChek Sweat Patch


- pH
- Conductivity
- Na⁺ ions concentration
- Temperature



A passive pump transports sweat to the multisensing patch



Acquisition and distribution layer (ADL)
to collect sweat from a larger surface





Integrated patch and electronics



Control Unit

pH Sensor

Patch Containing Sodium
and Conductivity Sensor



Jean Luprano, CSEM




Sodium Sensor



Na⁺ sensor: Ion Selective Electrode

Measurement of Open Circuit Potential (OCP) between Reference Gold Electrode and SC-ISE (Solid Contact- Ion Selective electrode)



With Dr. Isabelle Chartier CEA-LETI Grenoble


Na⁺ sensors (SC-ISE): patch testing

Calibration curves

NaCl model solutions $10 < [\text{NaCl}] < 100\text{mM}$

Linear variation OCP vs Log [Na⁺]
Good sensitivity 130-150mV/ log unit


Reversibility of Na⁺ ISE with model solutions



Ions	[ion] Hitachi 912 (mM)	[ion] ISE sensors (mM)	Typical concentration (literature, mM)
Na ⁺	53.8	50	54.5 ± 8.4

Measurements of sodium ions in
natural sweat
with two techniques (deviation 10%)


Na⁺ sensor: patch performance



Pictures of the experiment with a volunteer

(Experimental conditions : T = 45°C, hydrometry = 36%, speed 4km/h, slope 6%)

Real-time measurements



BUT

At the beginning, correlated measurements of OCP

Measured [Na⁺] ↗ when real [Na⁺] (Hitachi 912) ↘


- ⇒ Contamination of the pump ? Concentration effect in the pump ?
- ⇒ Evaporation or saturation of the pump ?

Na⁺ sensor: Sodium Sensor Belt (SSB)




Schazmann, B., et al., Analytical Methods, 2010, 2(4): p. 342-348.

SSB - CF trials






Literature
[Na⁺] Lower back
range: **26±19mM**



CF threshold > **60 mM** [Na⁺]

Na⁺ sensor: The Road Ahead





Na⁺ sensor: The Road Ahead





12.20 **A-2:L05** Recent Progress in Flexible Screen-printed Ion-selective Sensors for Environmental and Wearable Applications

G. MATZEU, G. ZULIANI, D. DIAMOND, CLARITY Centre for Sensor Web Technologies, NCSR, Dublin City University, Dublin, Ireland






pH Sensor




BIOTEX – pH sensor



Sweat pH: 5 - 7



Bromocresol Purple ($pK_a = 6.2$)



5.2

6.8

Emitter-detector LED's $\lambda = 660$ nm



Morris, D., et al., Sensors and Actuators B, 2009, **139**, p:231–236

BIOTEX – pH Experimental setup

Sensing patch and electronics incorporated onto a waist band

Reference



Fabric Patch

Hydrophobic surface


Waist band to secure sweat sensing system to reduce motion artefacts




Skincheck 1™ pH meter

BIOTEX – pH Trials Results

Subjects include healthy individuals and CF sufferers




(a) normal subject




(b) CF sufferer

Development of the Concept



Curto V., et al., Procedia Engineering, 2011, **25**, p:1561 – 1564




Main Characteristics:

- pH dye bromocresol purple
- smLED-photodiode detection
- Absorbent = passive pump

Curto V., et al., Sensors and Actuators B, 2012, doi: 10.1016/j.snb.2012.02.010

Three operative regimes




- **Dry state**
natural adsorption of the cotton thread
- **Capillarity**
absorbing capacity of the absorbent material
- **Saturation**
absorbent gets saturated

n=4

Average time to reach sensing area equal to 3 minutes

Micro-fluidic: Characterisation




	Basic t_d /min	Acidic t_d /min	Basic t_r /min	Acidic t_r /min
Cycle 1	4.74	0.60	4.96	1.7
Cycle 2	1.77	0.69	3.21	2.28
Cycle 3	2.37	1.13	3.97	3.06
Cycle 4	3.70	1.29	5.30	3.17

t_d – delay time

t_r – response time







Micro-fluidic: Performance



0 – 5 minutes warm-up



Lactate Sensor

Wearable Biosensing: New Trends




Prof. John Rogers @ University of Illinois

Organic ElectroChemical Transistors (OECTs)




PSS

PEDOT



$$V_g = 0$$



$$V_g \neq 0$$



Advantages



- ✓ Simple electrical readout
- ✓ Inherent signal amplification
- ✓ Easy incorporation into arrays and circuits
- ✓ Printing technologies make their fabrication particularly cost-effective for future industrialisation

Disadvantages


Need of an electrolyte to put in contact gate and channel

- ✓ complex liquid handling architecture
- ✓ leakage and contamination
- ✓ degradation of bio-receptor in solution (PBS)

Ionogel & OECTs: Lactate Sensor




- Hydrophilic IL - avoid
- Hydrated IL completely
- 20 μ L drop-casted on
- UV polymerisation:




lyte solution.

itation observed


Ionogel & OECTs: Lactate Sensor



25 mM Lactate




Compatible with detection of lactate in sweat (9 - 23 mM)



Lactate Sensor: toward wearability

- Novelty lies in the configuration of the sensor
- Solid State electrolyte on a flexible transistor based biosensor.
 - Implications for the wearability of the sensor



Conclusions

- Sensors based on textile greatly simplify the sampling of sweat and allow real-time measurements.
- Developments of ISE Na^+ sensor technology in order to enhance its wearability
- Improvements of real-time pH sensor
- Detection of lactate within relevant physiological range.
- Solid State electrolyte overcomes problems related with the design of the future wearable sensor

Acknowledgements



Prof. Dermot Dimond

Dr. Fernando Benito Lopez

Dr. Shirley Coyle

Giusy Matzeu

Dr. Claudio Zuliani



07/CE/I1147



Department of Bioelectronics
Microelectronics Center of Provence
www.emse.fr

Dion Khodagholy

Prof. George Malliaras

Prof. Roisin Owens



Thank you for your attention



Clinical trials with Cystic Fibrosis patients

- ❖ Cystic Fibrosis (CF) is the most common lethal genetic disease in Caucasians
- ❖ Exercise is also a useful method of increasing lung capacity and endurance in CF patients
- ❖ Low tolerance to climatic heat stress of sufferers might lead to increased morbidity and mortality
- ❖ pH has been reported to be more alkaline for subjects with CF; probably due to the increased sodium concentration
- ❖ pH/ Na^+ measurements might be useful for monitoring activity level to optimise physiotherapy / training regime